


NeuroPulse Master Clinical EMG Report

Patient: Patient Name | ID: NP-20260521-6079 | Date: 21-05-2026 16:41

1. Clinical Summary

Session Score	66/100
Symmetry Score	65.5%
Side Deficit	34.5% (Left lower)
Lowest Activity Channel	Channel 2
Fatigue Drop	Low (10.0% decline)
Signal Quality	Good

2. 30-Second Graph Snapshot



3. Muscle-wise Quantitative Summary

Muscle / Channel	Avg	Peak	Active %	Quality
Channel 1	0.9	2.0	0.0%	Good
Channel 2	0.5	0.7	0.0%	Low
Channel 3	0.8	1.2	0.0%	Good
Channel 4	1.3	2.1	0.0%	Good

4. Local Data-Based EMG Findings

■ Signal Reasoning (Quantitative Text Analysis):

- * Session Score: 67/100. Signal Quality: Good.
- * Side Deficit: 35.9%. Symmetry Score: 64.1%. Marked left-side deficit detected.
- * Lowest Activity: Channel 2. Highest Peak: Channel 4.
- * Fatigue Drop: Low (0.0% decline from early to recent activity).
- * Channel Balance: left-side channel balance 58%; right-side channel balance 49%.
- * Suggested focus: compare the two right-side muscles/channels and train the weaker activation pattern.
- * Selected affected side: Select Affected Side. Non-diagnostic rehab monitoring support only.

Recorded Session Findings:

- * Recording summary: 9790 EMG samples captured over approximately 60 seconds.
- * Session Score: 66/100. Signal Quality: Good.
- * Side Deficit: 34.5%. Symmetry Score: 65.5%. Impression: marked left-side activation deficit. Bilateral EMG symmetry is reduced and requires targeted rehab attention.
- * Lowest Activity: Channel 2 showed the lowest average activation (0.5).
- * Highest Peak: Channel 4 showed the highest peak activation (2.1).
- * Channel Balance: left-side channel balance 60.5%; right-side channel balance 58.5%.
- * Fatigue Drop: Low fatigue pattern with 10.0% decline from early to recent activity.
- * Sustained activation: Channel 1 0.0%, Channel 2 0.0%, Channel 3 0.0%, Channel 4 0.0% above threshold.
- * Clinical note: Affected side was not selected; interpretation is limited to channel activity and symmetry. This is a non-diagnostic rehab monitoring prototype that provides quantitative EMG feedback for physiotherapy support.

5. Clinical Observations & Notes

Summary Notes: N/A

Comparative Analysis:

- Left Side Avg: 0.7
- Right Side Avg: 1.0
- Symmetry Ratio (L/R): 65.5%
- Most Dynamic Channel: Channel 4